


YOUR
GLOBAL
PARTNER
IN
DIAGNOSTICS



virion\serion


SERION ELISA *classic* SERION ELISA *antigen*

Technische Daten / Technical Data


serion



multianalyt™




Artikelnr. Order No.	Artikel Product	Probe Sample	Verdünnung Dilution	Einheit Units	Messbereich Measurement Range	Grenzwertbereich Borderline Range	Packung Packaging
ESR 128 A	Adenovirus IgA	Serum / Plasma	1 + 500	U/ml	5 - 300	8 - 10 (Kind/Child); 11 - 14 (Erw./Adult)	96 Tests, 2 x STD, 1 CTRL
ESR 128 G	Adenovirus IgG	Serum / Plasma	1 + 2000	U/ml	1 - 180	8 - 10 (Kind/Child); 10 - 13 (Erw./Adult)	96 Tests, 2 x STD, 1 CTRL
ESR 132 A	Aspergillus fumigatus IgA	Serum / Plasma	1 + 100	U/ml	30 - 1000	50 - 70	96 Tests, 2 x STD, 1 CTRL
ESR 132 G	Aspergillus fumigatus IgG	Serum / Plasma	1 + 500	U/ml	30 - 1000	50 - 70	96 Tests, 2 x STD, 1 CTRL
ESR 132 M	Aspergillus fumigatus IgM	Serum / Plasma	1 + 500	U/ml	30 - 1000	50 - 70	96 Tests, 2 x STD, 1 CTRL
ESR 140 G	Bacillus anthracis IgG	Serum / Plasma	1 + 100	U/ml	5 - 600	10 - 15	96 Tests, 2 x STD, 1 CTRL
ESR 120 A	Bordetella pertussis IgA	Serum / Plasma	1 + 100	IU/ml	2 - 200	25 - 40	96 Tests, 2 x STD, 1 CTRL
ESR 120 G	Bordetella pertussis IgG	Serum / Plasma	1 + 100	IU/ml	10 - 1000	40 - 50	96 Tests, 2 x STD, 1 CTRL
ESR 120 M	Bordetella pertussis IgM	Serum / Plasma	1 + 100	U/ml	5 - 50	9 - 14	96 Tests, 2 x STD, 1 CTRL
ESR 1201 A	Bordetella pertussis Toxin IgA	Serum / Plasma	1 + 100	IU/ml	10 - 300	15 - 20	96 Tests, 2 x STD, 1 CTRL
ESR 1201 G	Bordetella pertussis Toxin IgG	Serum / Plasma	1 + 100	IU/ml	5 - 600	40 - 100	96 Tests, 2 x STD, 1 CTRL
ESR 121 G	Borrelia burgdorferi IgG	Serum / Plasma / Liquor (CSF)	1 + 100 ²⁾	U/ml	1 - 100	3 - 5	96 Tests, 2 x STD, 1 CTRL
ESR 121 M	Borrelia burgdorferi IgM	Serum / Plasma / Liquor (CSF)	1 + 100 ²⁾	U/ml	1 - 60	3 - 5	96 Tests, 2 x STD, 1 CTRL
ESR 116 A	Brucella IgA	Serum / Plasma	1 + 100	U/ml	5 - 100	10 - 15	96 Tests, 2 x STD, 1 CTRL
ESR 116 G	Brucella IgG	Serum / Plasma	1 + 100	U/ml	5 - 500	20 - 30	96 Tests, 2 x STD, 1 CTRL
ESR 116 M	Brucella IgM	Serum / Plasma	1 + 100	U/ml	5 - 100	15 - 20	96 Tests, 2 x STD, 1 CTRL
ESR 139 A	Campylobacter jejuni IgA	Serum / Plasma	1 + 1000	U/ml	5 - 1000	20 - 25	96 Tests, 2 x STD, 1 CTRL
ESR 139 G	Campylobacter jejuni IgG	Serum / Plasma	1 + 1000	U/ml	2 - 500	20 - 30	96 Tests, 2 x STD, 1 CTRL
ESR 139 M	Campylobacter jejuni IgM	Serum / Plasma	1 + 1000	U/ml	20 - 2000	40 - 60	96 Tests, 2 x STD, 1 CTRL
ESR 117 A	Candida IgA	Serum / Plasma	1 + 1000	U/ml	10 - 225	60 - 80	96 Tests, 2 x STD, 1 CTRL
ESR 117 G	Candida IgG	Serum / Plasma	1 + 1000	U/ml	5 - 450	40 - 100	96 Tests, 2 x STD, 1 CTRL
ESR 117 M	Candida IgM	Serum / Plasma	1 + 100	U/ml	5 - 500	60 - 80	96 Tests, 2 x STD, 1 CTRL
A 117	Candida	Serum	3 + 1	U/ml	0.8 - 57	1,15 - 2,3	96 Tests, 4 x STD, 2 CTRL
ESR 137 A	Chlamydia IgA	Serum / Plasma	1 + 100	U/ml	2 - 70	3 - 4	96 Tests, 2 x STD, 1 CTRL
ESR 137 G	Chlamydia IgG	Serum / Plasma	1 + 20	U/ml	5 - 800	10 - 15	96 Tests, 2 x STD, 1 CTRL
ESR 1371 A	Chlamydia pneumoniae IgA	Serum / Plasma	1 + 500	U/ml	4 - 150	10 - 13	96 Tests, 2 x STD, 1 CTRL
ESR 1371 G	Chlamydia pneumoniae IgG	Serum / Plasma	1 + 500	U/ml	4 - 200	10 - 12	96 Tests, 2 x STD, 1 CTRL
ESR 1371 M	Chlamydia pneumoniae IgM	Serum / Plasma	1 + 100	U/ml	3 - 200	10 - 15	96 Tests, 2 x STD, 1 CTRL
ESR 1372 A	Chlamydia trachomatis IgA	Serum / Plasma	1 + 100	U/ml	5 - 150	9 - 16	96 Tests, 2 x STD, 1 CTRL
ESR 1372 G	Chlamydia trachomatis IgG	Serum / Plasma	1 + 100	U/ml	4 - 250	10 - 15	96 Tests, 2 x STD, 1 CTRL
ESR 1372 M	Chlamydia trachomatis IgM	Serum / Plasma	1 + 100	U/ml	5 - 180	10 - 14	96 Tests, 2 x STD, 1 CTRL
ESR 1311 A	Coxiella burnetii Phase 1 IgA	Serum / Plasma	1 + 100		qualitativ / qualitative	OD C/O +/- 10 %	96 Tests, 1 x POS, 2 x C/O, 1 CTRL
ESR 1311 G	Coxiella burnetii Phase 1 IgG	Serum / Plasma	1 + 100		qualitativ / qualitative	OD C/O +/- 10 %	96 Tests, 1 x POS, 2 x C/O, 1 CTRL
ESR 1312 G	Coxiella burnetii Phase 2 IgG	Serum / Plasma	1 + 500	U/ml	5 - 500	20 - 30	96 Tests, 2 x STD, 1 CTRL
ESR 1312 M	Coxiella burnetii Phase 2 IgM	Serum / Plasma	1 + 100		qualitativ / qualitative	OD C/O +/- 10 %	96 Tests, 1 x POS, 2 x C/O, 1 CTRL



Artikelnr. Order No.	Artikel Product	Probe Sample	Verdünnung Dilution
ESR 134 A	Coxsackievirus IgA	Serum / Plasma	1 + 100
ESR 134 G	Coxsackievirus IgG	Serum / Plasma	1 + 500
ESR 134 M	Coxsackievirus IgM	Serum / Plasma	1 + 100
ESR 109 G	Cytomegalovirus IgG ^{A)}	Serum / Plasma / Liquor (CSF)	1 + 100
ESR 109 M	Cytomegalovirus IgM	Serum / Plasma	1 + 100
ESR 114 G	Dengue Virus IgG	Serum / Plasma	1 + 100
ESR 114 M	Dengue Virus IgM	Serum / Plasma	1 + 100
ESR 130 G	Diphtherie (Diphtheria) IgG	Serum / Plasma	1 + 100
ESR 107 G	Echinococcus IgG	Serum / Plasma	1 + 100
ESR 135 A	Echovirus IgA	Serum / Plasma	1 + 100
ESR 135 G	Echovirus IgG	Serum / Plasma	1 + 500
ESR 135 M	Echovirus IgM	Serum / Plasma	1 + 100
ESR 133 A	Enterovirus IgA	Serum / Plasma	1 + 100
ESR 133 G	Enterovirus IgG	Serum / Plasma	1 + 500
ESR 133 M	Enterovirus IgM	Serum / Plasma	1 + 100
ESR 1361 G	Epstein-Barr Virus VCA IgG	Serum / Plasma	1 + 100
ESR 1361 M	Epstein-Barr Virus VCA IgM	Serum / Plasma	1 + 100
ESR 1362 G	Epstein-Barr Virus EBNA1 IgG	Serum / Plasma	1 + 100
ESR 1363 G	Epstein-Barr Virus EA IgG	Serum / Plasma	1 + 100 ¹⁾
ESR 142 G	Francisella tularensis IgG	Serum / Plasma	1 + 100
ESR 142 M	Francisella tularensis IgM	Serum / Plasma	1 + 100
ESR 112 G	FSME (TBE) Virus IgG	Serum / Plasma / Liquor (CSF)	1 + 100
ESR 112 M	FSME (TBE) Virus IgM	Serum / Plasma / Liquor (CSF)	1 + 100
ESR 145 G	Hantavirus Puumala IgG	Serum / Plasma	1 + 100 ¹⁾
ESR 145 M	Hantavirus Puumala IgM	Serum / Plasma	1 + 100 ¹⁾
ESR 118 A	Helicobacter pylori IgA	Serum / Plasma	1 + 100
ESR 118 G	Helicobacter pylori IgG	Serum / Plasma	1 + 100
ESR 118 M	Helicobacter pylori IgM	Serum / Plasma	1 + 100
ESR 105 A	Herpes Simplex Virus 1+2 IgA	Serum / Plasma	1 + 100
ESR 105 G	Herpes Simplex Virus 1+2 IgG	Serum / Plasma / Liquor (CSF)	1 + 100
ESR 105 M	Herpes Simplex Virus 1+2 IgM	Serum / Plasma	1 + 100
ESR 1051 G	Herpes Simplex Virus 1 gG IgG	Serum / Plasma	1 + 100
ESR 1051 M	Herpes Simplex Virus 1 IgM	Serum / Plasma	1 + 100
ESR 1052 G	Herpes Simplex Virus 2 gG IgG	Serum / Plasma	1 + 100
ESR 1052 M	Herpes Simplex Virus 2 IgM	Serum / Plasma	1 + 100

Einheit Units	Messbereich Measurement Range	Grenzwertbereich Borderline Range	Packung Packaging
U/ml	5 - 500	30 - 50	96 Tests, 2 x STD, 1 CTRL
U/ml	15 - 1000	80 - 100	96 Tests, 2 x STD, 1 CTRL
U/ml	10 - 500	30 - 50	96 Tests, 2 x STD, 1 CTRL
PEI-U/ml	10 - 2000	25 - 40	96 Tests, 2 x STD, 1 CTRL
U/ml	5 - 600	10 - 15	96 Tests, 2 x STD, 1 CTRL
U/ml	5 - 500	10 - 15	96 Tests, 2 x STD, 1 CTRL
U/ml	5 - 200	10 - 15	96 Tests, 2 x STD, 1 CTRL
IU/ml	0,05 - 2,0		96 Tests, 2 x STD, 1 CTRL
U/ml	5 - 500	10 - 15	96 Tests, 2 x STD, 1 CTRL
U/ml	5 - 500	30 - 50	96 Tests, 2 x STD, 1 CTRL
U/m	15 - 1000	80 - 100	96 Tests, 2 x STD, 1 CTRL
U/ml	10 - 500	30 - 50	96 Tests, 2 x STD, 1 CTRL
U/ml	5 - 1000	30 - 50	96 Tests, 2 x STD, 1 CTRL
U/m	15 - 1000	80 - 100	96 Tests, 2 x STD, 1 CTRL
U/ml	5 - 1000	30 - 50	96 Tests, 2 x STD, 1 CTRL
U/ml	4 - 200	10 - 15	96 Tests, 2 x STD, 1 CTRL
U/ml	4 - 200	9 - 13	96 Tests, 2 x STD, 1 CTRL
U/ml	1 - 200	2,5 - 3,0	96 Tests, 2 x STD, 1 CTRL
U/ml	2 - 400	10 - 15	96 Tests, 2 x STD, 1 CTRL
U/ml	3 - 300	10 - 15	96 Tests, 2 x STD, 1 CTRL
U/ml	4 - 400	10 - 15	96 Tests, 2 x STD, 1 CTRL
U/ml	5 - 700	15 - 28	96 Tests, 2 x STD, 1 CTRL
U/ml	5 - 300	10 - 15	96 Tests, 2 x STD, 1 CTRL
U/ml	5 - 400	10 - 15	96 Tests, 2 x STD, 1 CTRL
U/ml	5 - 400	10 - 15	96 Tests, 2 x STD, 1 CTRL
U/ml	10 - 200	20 - 30	96 Tests, 2 x STD, 1 CTRL
U/ml	5 - 500	35 - 50	96 Tests, 2 x STD, 1 CTRL
U/ml	10 - 300	20 - 30	96 Tests, 2 x STD, 1 CTRL
U/ml	10 - 500	20 - 30	96 Tests, 2 x STD, 1 CTRL
U/ml	10 - 2000	20 - 30	96 Tests, 2 x STD, 1 CTRL
U/ml	10 - 500	20 - 30 / 60 - 70	96 Tests, 2 x STD, 1 CTRL
U/ml	10 - 1000	20 - 30	96 Tests, 2 x STD, 1 CTRL
U/ml	10 - 500	20 - 30	96 Tests, 2 x STD, 1 CTRL
U/ml	10 - 1000	20 - 30	96 Tests, 2 x STD, 1 CTRL
U/ml	10 - 500	20 - 30	96 Tests, 2 x STD, 1 CTRL



Artikelnr. Order No.	Artikel Product	Probe Sample	Verdünnung Dilution	Einheit Units	Messbereich Measurement Range	Grenzwertbereich Borderline Range	Packung Packaging
ESR 1231 A	Influenza A Virus IgA	Serum / Plasma	1 + 500	U/ml	2 - 200	6 - 9 (Kind/Child); 10 - 15 (Erw./Adult)	96 Tests, 2 x STD, 1 CTRL
ESR 1231 G	Influenza A Virus IgG	Serum / Plasma	1 + 2000	U/ml	2 - 200	4 - 7 (Kind/Child); 10 - 15 (Erw./Adult)	96 Tests, 2 x STD, 1 CTRL
ESR 1232 A	Influenza B Virus IgA	Serum / Plasma	1 + 500	U/ml	2 - 200	4 - 7 (Kind/Child); 10 - 15 (Erw./Adult)	96 Tests, 2 x STD, 1 CTRL
ESR 1232 G	Influenza B Virus IgG	Serum / Plasma	1 + 2000	U/ml	2 - 200	4 - 7 (Kind/Child); 10 - 15 (Erw./Adult)	96 Tests, 2 x STD, 1 CTRL
ESR 106 G	Legionella pneumophila 1-7 IgG	Serum / Plasma	1 + 100	U/ml	10 - 500	50 - 70	96 Tests, 2 x STD, 1 CTRL
ESR 106 M	Legionella pneumophila 1-7 IgM	Serum / Plasma	1 + 100	U/ml	10 - 500	120 - 140	96 Tests, 2 x STD, 1 CTRL
ESR 125 G	Leptospira IgG	Serum / Plasma	1 + 100	U/ml	2 - 100	10 - 15	96 Tests, 2 x STD, 1 CTRL
ESR 125 M	Leptospira IgM	Serum / Plasma	1 + 100	U/ml	3 - 100	15 - 20	96 Tests, 2 x STD, 1 CTRL
ESR 102 G	Masern (Measles) Virus IgG	Serum / Plasma / Liquor (CSF)	1 + 100	mIU/ml	50 - 5000	150 - 200	96 Tests, 2 x STD, 1 CTRL
ESR 102 M	Masern (Measles) Virus IgM	Serum / Plasma	1 + 100	U/ml	5 - 800	10 - 15	96 Tests, 2 x STD, 1 CTRL
ESR 103 G	Mumps (Parotitis) Virus IgG	Serum / Plasma / Liquor (CSF)	1 + 100	U/ml	30 - 2000	70 - 100	96 Tests, 2 x STD, 1 CTRL
ESR 103 M	Mumps (Parotitis) Virus IgM	Serum / Plasma	1 + 100	U/ml	5 - 600	10 - 15	96 Tests, 2 x STD, 1 CTRL
ESR 127 A	Mycoplasma pneumoniae IgA	Serum / Plasma	1 + 100	U/ml	2 - 150	10 - 14	96 Tests, 2 x STD, 1 CTRL
ESR 127 G	Mycoplasma pneumoniae IgG	Serum / Plasma	1 + 100	U/ml	3 - 200	10 - 15 (Kind/Child); 20 - 30 (Erw./Adult)	96 Tests, 2 x STD, 1 CTRL
ESR 127 M	Mycoplasma pneumoniae IgM	Serum / Plasma	1 + 100	U/ml	5 - 150	13 - 17	96 Tests, 2 x STD, 1 CTRL
ESR 126 A	Parainfluenza Virus 1, 2, 3 IgA	Serum / Plasma	1 + 100	U/ml	5 - 500	15 - 25	96 Tests, 2 x STD, 1 CTRL
ESR 1261 G	Parainfluenza Virus 1 IgG	Serum / Plasma	1 + 100	U/ml	5 - 600	25 - 35	96 Tests, 2 x STD, 1 CTRL
ESR 1262 G	Parainfluenza Virus 2 IgG	Serum / Plasma	1 + 100	U/ml	10 - 1000	90 - 150	96 Tests, 2 x STD, 1 CTRL
ESR 1263 G	Parainfluenza Virus 3 IgG	Serum / Plasma	1 + 100	U/ml	10 - 1000	65 - 95	96 Tests, 2 x STD, 1 CTRL
ESR 122 G	Parvovirus B19 IgG	Serum / Plasma	1 + 100	IU/ml	1,5 - 150	3 - 5	96 Tests, 2 x STD, 1 CTRL
ESR 122 M	Parvovirus B19 IgM	Serum / Plasma	1 + 100 ³⁾	U/ml	5 - 200	13 - 17	96 Tests, 2 x STD, 1 CTRL
ESR 113 A	Respiratory Syncytial Virus IgA	Serum / Plasma	1 + 500	U/ml	2 - 150	10 - 15	96 Tests, 2 x STD, 1 CTRL
ESR 113 G	Respiratory Syncytial Virus IgG	Serum / Plasma	1 + 2000	U/ml	2 - 200	10 - 15	96 Tests, 2 x STD, 1 CTRL
ESR 129 G	Röteln (Rubella) Virus IgG ^{A)}	Serum / Plasma / Liquor (CSF)	1 + 100	IU/ml	2 - 500	10 - 20	96 Tests, 2 x STD, 1 CTRL
ESR 129 M	Röteln (Rubella) Virus IgM	Serum / Plasma	1 + 100	U/ml	2 - 60	2,5 - 3,5	96 Tests, 2 x STD, 1 CTRL
ESR 108 G	Tetanus IgG	Serum / Plasma	1 + 100	IU/ml	0,05 - 5,0		96 Tests, 2 x STD, 1 CTRL
ESR 110 G	Toxoplasma gondii IgG ^{A)}	Serum / Plasma	1 + 100	IU/ml	5 - 500	10 - 20	96 Tests, 2 x STD, 1 CTRL
ESR 110 M	Toxoplasma gondii IgM	Serum / Plasma	1 + 100	U/ml	100 - 5000	300 - 350 / 450 - 540	96 Tests, 2 x STD, 1 CTRL
ESR 144 G	Thyroglobulin IgG	Serum / Plasma	1 + 100	IU/ml	20 - 2000	100 - 150	96 Tests, 2 x STD, 1 CTRL
ESR 143 G	Thyropoxidase IgG	Serum / Plasma	1 + 100	IU/ml	5 - 600	35 - 50	96 Tests, 2 x STD, 1 CTRL
ESR 104 A	Varicella-Zoster Virus IgA	Serum / Plasma	1 + 100	U/ml	5 - 1000	35 - 50	96 Tests, 2 x STD, 1 CTRL
ESR 104 G	Varicella-Zoster Virus IgG	Serum / Plasma / Liquor (CSF)	1 + 100	mIU/ml	15 - 2000	50 - 100	96 Tests, 2 x STD, 1 CTRL
ESR 104 M	Varicella-Zoster Virus IgM	Serum / Plasma	1 + 100	U/ml	5 - 1000	10 - 15	96 Tests, 2 x STD, 1 CTRL
ESR 138 A	Yersinia IgA	Serum / Plasma	1 + 100	U/ml	5 - 100	10 - 15	96 Tests, 2 x STD, 1 CTRL
ESR 138 G	Yersinia IgG	Serum / Plasma	1 : 20	U/ml	5 - 500	10 - 15	96 Tests, 2 x STD, 1 CTRL
ESR 138 M	Yersinia IgM	Serum / Plasma	1 + 100	U/ml	5 - 100	10 - 15	96 Tests, 2 x STD, 1 CTRL

Unsere SERION ELISA *classic* und SERION ELISA *antigen* sind validiert für den manuellen Gebrauch sowie für die Anwendung auf Immunomat™, DYNEX DSX™ und DS2.

SERION ELISA *classic* and SERION ELISA *antigen* are evaluated for manual testing as well as on the Immunomat™, DYNEX DSX™ and DS2.



Institut Virion\Serion GmbH

Tel./Phone +49 931 3045-0
Fax +49 931 3045-100

Serion Immundiagnostica GmbH

Tel./Phone +49 931 3045-580
Fax +49 931 3045-590

Friedrich-Bergius-Ring 19, D-97076 Würzburg, Germany
Postfach 92 53, D-97076 Würzburg, Germany
E-Mail: dialog@virion-serion.de Internet: www.virion-serion.de